

Anti-PTPMT1 antibody (1-201) (STJ111326)

STJ111326

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PTPMT1 (1-201) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	114971
Gene Symbol	PTPMT1
Uniprot ID	PTPM1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-201 of human PTPMT1 (NP_783859.1).
Immunogen Region	1-201
Specificity	
Immunogen Sequence	