

Anti-A4GALT antibody (244-353) (STJ111431)

STJ111431

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-A4GALT (244-353) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53947
Gene Symbol	A4GALT
Uniprot ID	A4GAT_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 244-353 of human A4GALT (NP_059132.1).
Immunogen Region	244-353
Specificity	
Immunogen Sequence	