

Anti-SCAMP1 antibody (1-150) (STJ111561)

STJ111561

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-SCAMP1 (1-150) is suitable for use in Western Blot.

Applications WB
Host/Source Rabbit

Reactivity Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID

Concentration

Conjugation
Purification
Dilution Range
Unconjugated
Affinity purification
WB 1:500-1:2000

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

Isotype IgG


Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID 9522
Gene Symbol SCAMP1
Uniprot ID SCAM1_HUMAN

Immunogen Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human SCAMP1 (NP_004857.4).

Immunogen Region 1-150 Specificity Immunogen Sequence



Western blot analysis of extracts of various cell line using SCAMP1 antibody (STJ111561) at 1:100 dilution. Secondary antibody: HRP Goat Anti-rabbit Ig (H+L) at 1:10000 dilution. Lysates/proteins: 25ug price lane. Blocking buffer: 3% nonfat dry milk in TBS