

Anti-PDGFR α antibody (15-205) (STJ111596)

STJ111596

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PDGFR α (15-205) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5157
Gene Symbol	PDGFRα
Uniprot ID	PDGFR_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 15-205 of human PDGFR α (NP_006198.1).
Immunogen Region	15-205
Specificity	
Immunogen Sequence	