

Anti-WRNIP1 antibody (401-640) (STJ111597)

STJ111597

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-WRNIP1 (401-640) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56897
Gene Symbol	WRNIP1
Uniprot ID	WRIP1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 401-640 of human WRNIP1 (NP_569079.1).
Immunogen Region	401-640
Specificity	
Immunogen Sequence	