

Anti-COL9A3 antibody (445-684) (STJ111840)

STJ111840

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-COL9A3 (445-684) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:100-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1299
Gene Symbol	COL9A3
Uniprot ID	C09A3_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 445-684 of human COL9A3 (NP_001844.3).
Immunogen Region	445-684
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded mouse heart using COL9A3 antibody (STJ111840) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081