

Anti-KCNJ9 antibody (1-100) (STJ111867)

STJ111867

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-KCNJ9 (1-100) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3765
Gene Symbol	KCNJ9
Uniprot ID	KCNJ9_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KCNJ9 (NP_004974.2).
Immunogen Region	1-100
Specificity	
Immunogen Sequence	