

Anti-MYBL1 antibody (443-692) (STJ111871)

STJ111871

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MYBL1 (443-692) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4603
Gene Symbol	MYBL1
Uniprot ID	MYBA_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 443-692 of human MYBL1 (NP_001138227.1).
Immunogen Region	443-692
Specificity	
Immunogen Sequence	