

Anti-NFYC antibody (86-335) (STJ111874)

STJ111874

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NFYC (86-335) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4802
Gene Symbol	NFYC
Uniprot ID	NFYC_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 86-335 of human NFYC (NP_055038.2).
Immunogen Region	86-335
Specificity	
Immunogen Sequence	