

Anti-DSTYK antibody (635-884) (STJ111912)

STJ111912

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DSTYK (635-884) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25778
Gene Symbol	DSTYK
Uniprot ID	DUSTY_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 635-884 of human DSTYK (NP_955749.1).
Immunogen Region	635-884
Specificity	
Immunogen Sequence	