

Anti-RPP25L antibody (1-163) (STJ111931)

STJ111931

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RPP25L (1-163) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	138716
Gene Symbol	RPP25L
Uniprot ID	RP25L_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human RPP25L (NP_680544.1).
Immunogen Region	1-163
Specificity	
Immunogen Sequence	