

Anti-EV71 antibody (STJ112043)

STJ112043

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-EV71 3C is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Enterovirus

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:5000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant protein of Enterovirus EV71 3C
Immunogen Region	
Specificity	
Immunogen Sequence	