

Anti-TNFSF9 antibody (STJ112240)

STJ112240

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-4-1BBL/CD137L is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8744
Gene Symbol	TNFSF9
Uniprot ID	TNFL9_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-254 of human 4-1BBL/CD137L (NP_003802.1).
Immunogen Region	
Specificity	
Immunogen Sequence	