

Anti-TST antibody (1-297) (STJ112560)

STJ112560

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TST (1-297) is suitable for use in Western Blot, Immunohistochemistry and Immunoprecipitation.
Applications	WB, IHC, IP
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000 IHC 1:50-1:200 IP 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7263
Gene Symbol	TST
Uniprot ID	THTR_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-297 of human TST (NP_003303.2).
Immunogen Region	1-297
Specificity	
Immunogen Sequence	




Western blot analysis of extracts of various cell lines, using TST antibody (STJ112560) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.



Immunohistochemistry of paraffin-embedded rat kidney using TST rabbit polyclonal antibody (STJ112560) at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded rat testis using TST rabbit polyclonal antibody (STJ112560) at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human liver using TST rabbit polyclonal antibody (STJ112560) at dilution of 1:200 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081