

Anti-SEC14L2 antibody (60-270) (STJ112600)

STJ112600

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SEC14L2 (60-270) is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence.
Applications	WB, IHC, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000 IHC 1:50-1:200 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION


Gene ID	23541
Gene Symbol	SEC14L2
Uniprot ID	S14L2_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-270 of human SEC14L2 (NP_001191133.1).
Immunogen Region	60-270
Specificity	
Immunogen Sequence	




Western blot analysis of extracts of various cell lines, using SEC14L2 antibody (STJ112600) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunohistochemistry of paraffin-embedded human esophageal cancer using SEC14L2 rabbit polyclonal antibody (STJ112600) at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded mouse pancreas using SEC14L2 rabbit polyclonal antibody (STJ112600) at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded rat lung using SEC14L2 rabbit polyclonal antibody (STJ112600) at dilution of 1:200 (40x lens).