

Anti-NARS antibody (309-548) (STJ113269)

STJ113269

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NARS (309-548) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4677
Gene Symbol	NARS1
Uniprot ID	SYNC_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 309-548 of human NARS (NP_004530.1).
Immunogen Region	309-548
Specificity	
Immunogen Sequence	