

Anti-GABRA3 antibody (29-276) (STJ113274)

STJ113274

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GABRA3 (29-276) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2556
Gene Symbol	GABRA3
Uniprot ID	GBRA3_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-276 of human GABRA3 (NP_000799.1).
Immunogen Region	29-276
Specificity	
Immunogen Sequence	