

Anti-ATG7 antibody (500-676) (STJ113569)

STJ113569

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ATG7 (500-676) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10533
Gene Symbol	ATG7
Uniprot ID	ATG7_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-676 of human ATG7 (NP_001129503.2).
Immunogen Region	500-676
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded human breast cancer using ATG7 antibody (STJ113569) at dilution of 1:100 (40x lens).