

Anti-C1QTNF4 antibody (250-329) (STJ113603)

STJ113603

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-C1QTNF4 (250-329) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	114900
Gene Symbol	C1QTNF4
Uniprot ID	C1QT4_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-329 of human C1QTNF4 (NP_114115.2).
Immunogen Region	250-329
Specificity	
Immunogen Sequence	