

Anti-IL1B antibody (124-223) (STJ113778)

STJ113778

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IL1 beta (124-223) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3553
Gene Symbol	IL1B
Uniprot ID	IL1B_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 124-223 of human IL1 Beta (NP_000567.1).
Immunogen Region	124-223
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of THP-1 cells using IL1B antibody (STJ113778) at dilution of 1:100. Blue: DAPI for nuclear staining.