

Anti-SNRPA1 antibody (1-220) (STJ114054)

STJ114054

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SNRPA1 (1-220) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6627
Gene Symbol	SNRPA1
Uniprot ID	RU2A_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human SNRPA1 (NP_003081.2).
Immunogen Region	1-220
Specificity	
Immunogen Sequence	




Western blot analysis of extracts of various cell lines, using SNRPA1 antibody (STJ114054) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunohistochemistry of paraffin-embedded human esophageal cancer using SNRPA1 rabbit polyclonal antibody (STJ114054) at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human lung cancer using SNRPA1 rabbit polyclonal antibody (STJ114054) at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded mouse liver using SNRPA1 rabbit polyclonal antibody (STJ114054) at dilution of 1:200 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081