

## Anti-PCSK1 antibody (604-753) (STJ114128) STJ114128

## **GENERAL INFORMATION**

 Product Type
 Primary antibodies

 Short Description
 Rabbit polyclonal antibody anti-PCSK1 (604-753) is suitable for use in Western Blot and Immunofluorescence.

 Applications
 WB, IF

 Host/Source
 Rabbit

 Reactivity
 Mouse, Rat

## **PRODUCT PROPERTIES**

 
 Clonality Clone ID
 Polyclonal

 Concentration

 Conjugation
 Uconjugated


 Purification
 Affinity purification

 Dilution Range
 WB 1:500-1:2000 IF 1:50-1:200

 Formulation
 PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3. Isotype

 IgG
 Store in a freezer at-20°C and avoid freeze-thaw cycles.

## **TARGET INFORMATION**



Western blot analysis of extracts of rat brain, using PCSK1 antibody (STJ114128) at 1:3000 dilution. Secondary antibody. HPC Goat Anti-rabbit IGG (H+1) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic KT. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes. St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081