

Anti-ELN antibody (STJ114307)

STJ114307

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ELN is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2006
Gene Symbol	ELN
Uniprot ID	ELN_HUMAN
Immunogen	Recombinant protein of human ELN
Immunogen Region	
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded human liver cancer using ELN antibody (STJ114307) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human breast using ELN antibody (STJ114307) at dilution of 1:100 (40x lens).