

Anti-SUMO1 antibody (STJ114403)

STJ114403

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SUMO1 is suitable for use in Immunohistochemistry and Immunofluorescence.
Applications	IHC, IF
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7341
Gene Symbol	SUMO1
Uniprot ID	SUMO1_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 20 to the C-terminus of human SUMO1 (NP_003343.1).
Immunogen Region	
Specificity	
Immunogen Sequence	