

Anti-IKBKG antibody (40-140) (STJ114435)

STJ114435

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IKK Gamma (40-140) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8517
Gene Symbol	IKBKG
Uniprot ID	NEMO_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 40-140 of human IKK Gamma (NP_003630.1).
Immunogen Region	40-140
Specificity	
Immunogen Sequence	