

Anti-MOAP1 antibody (1-351) (STJ114473)

STJ114473

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MOAP1 (1-351) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64112
Gene Symbol	MOAP1
Uniprot ID	MOAP1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-351 of human MOAP1 (NP_071434.2).
Immunogen Region	1-351
Specificity	
Immunogen Sequence	