

Anti-AGPAT5 antibody (120-320) (STJ114750)

STJ114750

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-AGPAT5 (120-320) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55326
Gene Symbol	AGPAT5
Uniprot ID	PLCE_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-320 of human AGPAT5 (NP_060831.2).
Immunogen Region	120-320
Specificity	
Immunogen Sequence	