

Anti-BMF antibody (1-184) (STJ114773)

STJ114773

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-BMF (1-184) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	90427
Gene Symbol	BMF
Uniprot ID	BMF_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-184 of human BMF (NP_277038.1).
Immunogen Region	1-184
Specificity	
Immunogen Sequence	