

Anti-Phospho-ARAF-S299 antibody (STJ114847)

STJ114847

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-ARAF-S299 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	369
Gene Symbol	ARAF
Uniprot ID	ARAF_HUMAN
Immunogen	A synthetic phosphorylated peptide around S299 of human ARAF (NP_001645.1).
Immunogen Region	
Specificity	
Immunogen Sequence	