

Anti-USP6 antibody (55-120) (STJ115220)

STJ115220

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-USP6 (55-120) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9098
Gene Symbol	USP6
Uniprot ID	UBP6_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-120 of human USP6 (NP_004496.2).
Immunogen Region	55-120
Specificity	
Immunogen Sequence	