

Anti-MAP4K3 antibody (755-894) (STJ115337)

STJ115337

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-MAP4K3 (755-894) is suitable for use in Immunofluorescence.

Applications IF
Host/Source Rabbit
Reactivity Human

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID Concentration

Conjugation
Purification
Dilution Range
Unconjugated
Affinity purification
F1:50-1:200

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

Isotype IgG

Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID 8491 Gene Symbol MAP4K3 Uniprot ID M4K3_HUMAN

Immunogen Recombinant fusion protein containing a sequence corresponding to amino acids 755-894 of human MAP4K3

(NP_003609.2).

Immunogen Region 755-894

Specificity Immunogen Sequence