

Anti-PIAS3 antibody (399-628) (STJ115356)

STJ115356

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PIAS3 (399-628) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10401
Gene Symbol	PIAS3
Uniprot ID	PIAS3_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 399-628 of human PIAS3 (NP_006090.2).
Immunogen Region	399-628
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of U2OS cells using PIAS3 antibody (STJ115356).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081