

Anti-UBD antibody (1-165) (STJ115359)

STJ115359

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-UBD (1-165) is suitable for use in Immunohistochemistry and Immunofluorescence.
Applications	IHC, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10537
Gene Symbol	UBD
Uniprot ID	UBD_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human FAT10/FAT10/UBD (NP_006389.2).
Immunogen Region	1-165
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of A549 cells using FAT10/FAT10/UBD antibody (STJ115359). Blue: DAPI for nuclear staining.