

Anti-SETD6 antibody (145-410) (STJ115400)

STJ115400

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SETD6 (145-410) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79918
Gene Symbol	SETD6
Uniprot ID	SETD6_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 145-410 of human SETD6 (NP_001153777.1).
Immunogen Region	145-410
Specificity	
Immunogen Sequence	