

Anti-BCL2L2 antibody (1-193) (STJ115432)

STJ115432

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Bcl-W (1-193) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	599
Gene Symbol	BCL2L2
Uniprot ID	B2CL2_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human Bcl-W (NP_001186768.1).
Immunogen Region	1-193
Specificity	
Immunogen Sequence	