

Anti-P4HTM antibody (294-563) (STJ115554)

STJ115554

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-P4HTM (294-563) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54681
Gene Symbol	P4HTM
Uniprot ID	P4HTM_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 294-563 of human P4HTM (NP_808807.2).
Immunogen Region	294-563
Specificity	
Immunogen Sequence	