

Anti-MSMB antibody (21-114) (STJ115584)

STJ115584

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MSMB (21-114) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4477
Gene Symbol	MSMB
Uniprot ID	MSMB_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-114 of human MSMB (NP_002434.1).
Immunogen Region	21-114
Specificity	
Immunogen Sequence	