

Anti-CELA3B antibody (80-190) (STJ115668)

STJ115668

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CELA3B (80-190) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23436
Gene Symbol	CELA3B
Uniprot ID	CEL3B_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-190 of human CELA3B (NP_031378.1).
Immunogen Region	80-190
Specificity	
Immunogen Sequence	



Western blot analysis of extracts of rat pancreas, using CELA3B antibody (STJ115668) at 1:3000 dilution.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081