

Anti-ADCYAP1R1 antibody (21-155) (STJ115798)

STJ115798

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ADCYAP1R1 (21-155) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	117
Gene Symbol	ADCYAP1R1
Uniprot ID	PACR_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-155 of human ADCYAP1R1 (NP_001186564.1).
Immunogen Region	21-155
Specificity	
Immunogen Sequence	