

Anti-CYP1B1 antibody (244-543) (STJ115902)

STJ115902

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CYP1B1 (244-543) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1545
Gene Symbol	CYP1B1
Uniprot ID	CP1B1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 244-543 of human CYP1B1 (NP_000095.2).
Immunogen Region	244-543
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of HeLa cells using CYP1B1 antibody (STJ115902) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.