

Anti-EDN3 antibody (17-238) (STJ115909)

STJ115909

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-EDN3 (17-238) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1908
Gene Symbol	EDN3
Uniprot ID	EDN3_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-238 of human EDN3 (NP_996917.1).
Immunogen Region	17-238
Specificity	
Immunogen Sequence	