

Anti-PLOD2 antibody (470-650) (STJ115980)

STJ115980

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PLOD2 (470-650) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5352
Gene Symbol	PLOD2
Uniprot ID	PLOD2_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 470-650 of human PLOD2 (NP_891988.1).
Immunogen Region	470-650
Specificity	
Immunogen Sequence	