

Anti-RPA3 antibody (6-121) (STJ115993)

ST.1115993

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-RPA3 (6-121) is suitable for use in Western Blot and Immunofluorescence.

Applications WB, IF Host/Source Rabbit

Reactivity Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID Concentration

Conjugation
Purification Unconjugated
Affinity purification
Dilution Range WB 1:500-1:2000

IF 1:50-1:200

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

Isotype IgG

Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.


TARGET INFORMATION

Gene ID 6119
Gene Symbol RPA3


Uniprot ID RFA3_HUMAN

Immunogen Recombinant fusion protein containing a sequence corresponding to amino acids 6-121 of human RPA3 (NP_002938.1).

Immunogen Region 6-121 Specificity Immunogen Sequence



Western blot analysis of extracts of mouse spleen, using PRA3 antibody (STJ115993) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 39% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exoosure time: 60s.



Immunofluorescence analysis of MCF7 cells using RPA3 antibody (STJ115993).