

Anti-HYAL3 antibody (198-417) (STJ116010)

STJ116010

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HYAL3 (198-417) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8372
Gene Symbol	HYAL3
Uniprot ID	HYAL3_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 198-417 of human HYAL3 (NP_003540.2).
Immunogen Region	198-417
Specificity	
Immunogen Sequence	