

Anti-RNF8 antibody (1-260) (STJ116018)

ST.1116018

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-RNF8 (1-260) is suitable for use in Western Blot and Immunofluorescence.

Applications WB, IF Host/Source Rabbit Reactivity Human

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID Concentration

Conjugation
Purification
Dilution Range
Unconjugated
Affinity purification
WB 1:500-1:2000

IF 1:50-1:200

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

Isotype IgG

Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID 9025
Gene Symbol RNF8


Uniprot ID RNF8_HUMAN

Immunogen Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human RNF8 (NP_003949.1).

Immunogen Region 1-260

Specificity

Immunogen Sequence



Immunofluorescence analysis of MCF-7 cells using RNF8 antibody (STJ116018). Blue: DAPI for nuclear