

Anti-KCNE2 antibody (1-123) (STJ116033)

STJ116033

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-KCNE2 (1-123) is suitable for use in Western Blot. |
| Applications | WB |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | WB 1:500-1:2000 |
| Formulation | PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3. |
| Isotype | IgG |
| Storage Instruction | Store in a freezer at -20°C and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 9992 |
| Gene Symbol | KCNE2 |
| Uniprot ID | KCNE2_HUMAN |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 1-123 of human KCNE2 (NP_751951.1). |
| Immunogen Region | 1-123 |
| Specificity | |
| Immunogen Sequence | |