

## Anti-EHF antibody (1-210) (STJ116050)

STJ116050

## **GENERAL INFORMATION**

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-EHF (1-210) is suitable for use in Western Blot.

Applications WB
Host/Source Rabbit
Reactivity Mouse

## **PRODUCT PROPERTIES**

Clonality Polyclonal

Clone ID

Concentration

Conjugation
Purification
Dilution Range
Unconjugated
Affinity purification
WB 1:500-1:2000

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

**Isotype** IgG

Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.

## **TARGET INFORMATION**


Gene ID 26298 Gene Symbol EHF

Uniprot ID EHF\_HUMAN

Immunogen Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human EHF (NP\_036285.2).

Immunogen Region 1-210

Specificity Immunogen Sequence



Western blot analysis of extracts of mouse lung, using EHF antibody (STJ116050) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer 3% honfat dry milk in TBST. Detection: