

## Anti-GAS7 antibody (1-412) (STJ116112)

## **GENERAL INFORMATION**

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-GAS7 (1-412) is suitable for use in Western Blot and Immunofluorescence.

Applications WB, IF Host/Source Rabbit

Reactivity Human, Mouse, Rat

## **PRODUCT PROPERTIES**

Clonality Polyclonal

Clone ID Concentration

Conjugation Unconjugated Purification Affinity purification Dilution Range WB 1:500-1:2000

IF 1:50-1:200

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

Isotype IgG

Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.

## **TARGET INFORMATION**

Gene ID 8522 Gene Symbol GAS7 Uniprot ID GAS7\_HUMAN

Immunogen Recombinant fusion protein containing a sequence corresponding to amino acids 1-412 of human GAS7 (NP\_001124303.1).

Immunogen Region 1-412 Specificity Immunogen Sequence



