

Anti-TAF1B antibody (369-588) (STJ116134)

STJ116134

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TAF1B (369-588) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9014
Gene Symbol	TAF1B
Uniprot ID	TAF1B_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 369-588 of human TAF1B (NP_005671.2).
Immunogen Region	369-588
Specificity	
Immunogen Sequence	