

Anti-CYP3A4 antibody (244-503) (STJ116145)

STJ116145

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CYP3A4 (244-503) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1576
Gene Symbol	CYP3A4
Uniprot ID	CP3A4_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 244-503 of human CYP3A4 (NP_059488.2).
Immunogen Region	244-503
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of HeLa cells using CYP3A4 antibody (STJ116145) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.